



www.epa.gov

Characterization of Nanoparticles at Environmentally Relevant Concentrations in Aquatic Media by Coupling Asymmetric Flow Field-Flow Fractionation with Single Particle – Inductively Coupled Plasma Mass Spectrometry

[†]National Research Council Post-Doctoral Associate at U.S. EPA, National Exposure Research Laboratory, Environmental Sciences Division, Las Vegas, NV 89119

[‡]Chemistry Department, University of Nevada Las Vegas, NV 89154

[§]U.S. EPA, National Exposure Research Laboratory, Environmental Sciences Division, Las Vegas, NV 89119

Khanh An Huynh,^{†*} Emily Siska,[‡] and Edward Heithmar[§]

Abstract

An advanced analytical method using asymmetric flow field-flow fractionation (AsFIFFF) and inductively coupled plasma mass spectrophotometry operated in single particle mode (spICP-MS) is developed to detect and quantify silver nanoparticles (AgNPs). Different from previous studies, the AsFIFFF and spICP-MS are directly connected for real-time determination of mass-based and hydrodynamic diameters of nanoparticles. Experimental results show that the AsFIFFF-spICP-MS is capable of detecting and quantifying AgNPs in complex mixtures such as the one containing 40, 60, 80 nm AgNPs and Ag-SiO₂ core-shell nanoparticles

Objectives

- Optimizing the performance of AsFIFFF in separating 40 – 100 nm AgNPs within 30 min: carrier, channel flow, cross flow, focusing time, particle concentrations, and cleaning procedures
- Optimizing spICP-MS measurements: dwell-time and particle concentrations
- Investigate the performance of the AsFIFFF-spICP-MS using mixtures containing different sizes and kinds of AgNPs

Background

- AgNPs are the most commonly used engineered nanomaterials in consumer products → AgNPs are likely to end up in the environment and possibly cause adverse impacts on ecosystem and human health
- To study the fate and transport, and assess the risk of AgNPs, it is critical to detect and quantify AgNPs in environmental samples
- Currently, methods used to detect, quantify, and characterize of AgNPs are not well developed
- AsFIFFF (combining with ICP-MS) and spICP-MS alone cannot provide detailed AgNP information such as coating layer thickness and aggregation state
- By using AsFIFFF-spICP-MS, the coating layer thickness and aggregation state of AgNPs can be determined since AsFIFFF-spICP-MS is able to measure both hydrodynamic diameter and mass-based diameter of detected AgNPs simultaneously.




www.silverlon.com/adhesive_strip.html



store.sealshield.com/silver-seal-keyboard-p26.aspx

Coupling AsFIFFF with sp-ICP-MS



The channel flow from the AsFIFFF directly entered the ICP-MS

AF2000 Focus (Postnova Analytic)

- Carrier 0.02 % FL-70
- Membrane 10 kDa RC
- Channel thickness 350 μ m
- Channel flow 1.5 mL/min
- Cross flow 1.2 mL/min
- Focusing time 3 min


NexION 300D (PerkinElmer) with Nano Application for Syngistix ICPMS Software

- Flow rate 1.5 mL/min
- Dwell-time 5 ms
- Sampling time 30 min
- Data acquisition after focusing


Nanoparticles

- 60 nm AuNPs (NIST 8013)
- 40, 60, 80, and 100 nm AgNPs (nanoComposix)
- Ag-SiO₂ core-shell nanoparticles (nanoComposix, total diameter = 92.3 nm, Ag core diameter = 51.0 nm) – modeled (i) AgNPs with coating layer or (ii) AgNP heteroaggregates



60 nm AgNPs, TEM from nanoComposix



Ag-SiO₂ NPs, TEM from nanoComposix






Calibration curves used in AsFIFFF-spICP-MS




AsFIFFF-spICP-MS Results


60 nm AgNPs and Ag-SiO₂ NPs



40 – 80 nm AgNPs and Ag-SiO₂ NPs



Recovery and Reproducibility



Conclusions. AsFIFFF-spICP-MS successfully (i) differentiates 60 nm AgNPs from Ag-SiO₂ core-shell NPs and (ii) detects and quantifies AgNPs in a complex mixture containing 40 – 60 nm AgNPs and Ag-SiO₂ core-shell NPs. The recovery and reproducibility of this technique are good to moderate depending on particle size and type.

Acknowledgement. We would like to thank Dr. Soheyl Tadjiki at Postnova for FFF training

Notice: the U.S. EPA, through its Office of Research and Development (ORD) funded this research and approved this abstract as a basic for a poster presentation. The actual presentation has not been peer reviewed by EPA. Mention of trade names or commercial products does not constitute endorsement or recommendation by EPA for use

***Email:** huynh.khanhan@epa.gov