

# Land Application Field Study II



**EPA: Ron Herrmann, Carolyn Acheson, Larry Zintek, Danielle Kleinmaier, Amanda Wroble, Greg Mitsakopoulos, and Kathleen Swan**

**PTSI: Andrea Burkes, Josh Kickish, Bob Grosser, Babina Shrestha**

# Outline

- Previous Research
- Experimental Design
- Data
  - Metals
  - Nonylphenol
  - PFAS
- Conclusions
- Next Steps
- Acknowledgements





# Disclaimer

**The views expressed in this presentation are those of the authors and do not necessarily reflect the views or policies of the U.S. Environmental Protection Agency.**



# Multimedia Land Application Study: Field Study I

- Surface application by side discharge manure spreader
- Agronomic rate of 10 wet tons/acre
- Material applied
  - Anaerobically digested biosolids
  - Polymer addition during dewatering
  - Lime addition
- Application field
  - Fescue field
  - No prior application of biosolids
  - Autumn application
  - Sampled for 1 month before and 4 months after application






# Soil Study Activities

- Characterize Study Conditions
  - Weather data
  - Soil data
  - Quantity and distribution of biosolids
  - Microbial community quantity and structure
- Performance Measurements
  - Microbes: fecal coliform density, viable helminth ova, Salmonella, enteric viruses, coliphage
  - Chemicals: concentrations of alkylphenol ethoxylates and degradation products (APEs)
  - Ecotoxicity Screening




## Soil Study Conclusions

- Changes observed in shallow samples after application
  - Microbial community shifted for about 28 days after application
  - Total biomass, fecal coliforms, and APEs
    - Increased following application
    - Persisted for 98 day sampling period
- 
- See full results in report “Multimedia Sampling During the Application of Biosolids on a Land Test Site”
    - [Report - https://www.epa.gov/sites/production/files/2018-11/documents/multimedia-sampling-land-testsite.pdf](https://www.epa.gov/sites/production/files/2018-11/documents/multimedia-sampling-land-testsite.pdf)
    - Summary - <https://www.epa.gov/sites/production/files/2018-11/documents/study-examines-fate-agricultural-land.pdf>

# Land Application Field Study II

- Research Questions
  - How does concentration change with time when biosolids are land applied?
  - Does the application method (Solid or Liquid) affect measured concentrations?
- Pilot/Field scale Treatment plot at local WWTP on a fescue and rye grass field
- Fall application at 10 wet tons/acre
- Study Design
  - Land application techniques (liquid and solids)\*
  - No application (control) and biosolids only (blue circles)\*
  - 3 Treatment reps of each
  - Sampled for 13 months\*
- Analytes
  - Microbes: fecal coliforms, total biomass and community structure
  - Nutrients
  - Chemicals: metals, APEs, and PFAS\*





## Samples from the Plots



**Control**

**Solid**


**Liquid**



## Sodium Data

### Concentrations


- Elevated in the solids and liquid treatments after application
- By day 120 near control levels



# Copper Data


## Concentrations

- Higher in the solid treatment throughout the study
- Liquid and control similar



# Nonylphenol (NP) Data

- Aerobically degradable surfactant, weakly estrogenic
- Only concentrations above the reporting limit (RL) are shown
- RL ~ 120  $\mu\text{g/kg}$  dry mass
- Liquid – no data > RL after 120 days
- NP persists in solid and biosolids throughout the study
- 34 % of samples did not meet QA acceptance criteria






# PFAS by ASTM D7968 (LC/MS/MS )

- Matrix - Environmental solids such as soils, sediments, and sludges
  - Developed by Larry Zintek (Reg 5 Chicago Regional Laboratory)
  - Single lab validated
- Method
  - Solvent extraction
  - Analysis by LC/MS/MS with MRMs and ion ratios
- Target Analytes:
  - 11 Perfluorinated Carboxylic Acids (PFCAs): C4 - C14
  - 3 Perfluorinated Sulfonic Acids (PFSAs): C4-C10
  - Precursors
    - 6 PFCAs - 6:2, 8:2, 10:2, & 7:3 FTCA; and 6:2 & 8:2 FTUCA
- Surrogate standards (isotopically labeled compounds): 9 PFCAs and PFSAs
  - Used to monitor analytical method performance/quality
  - Not used to “correct” the data




# LC/MS/MS Analytical Method – ASTM D7968



# Analytical Method Quality Controls


- Analyte Identification
  - Each batch: Initial calibration, Calibration check, and Second source check
  - Each analyte: Retention time, Primary and Confirmation ion masses, and Ion ratio
- Accuracy – 2 of each per batch unless specified
  - Surrogate spiking - All samples and blanks
    - Used to assess method performance
    - Not used to alter reported concentrations
  - Matrix spike samples – MS and MS duplicates
  - Spiked blanks
  - Method reporting limit checks
- Precision - 2 of each per batch
  - Duplicate samples
  - Matrix spike duplicates
  - Spiked blanks
- Laboratory Contamination – method blanks – 2 per batch





# ASTM D7968 Performance Data

**Mean  
Recovery  
(%)**




- Error bars are % RSD
- 6 replicates of each matrix
- Spiked at 400 ng/kg dry soil for all except 8:2 FTCA 8000 ng/kg dry soils
- 4 ASTM soil matrices: CL-1; CH-1; SP-1; ML-16

\* PFOS not shown for SP-1 and ML-1 because the matrices had background conc comparable to spike conc


# PFAS with EPA Screening Levels

- Conc above RL are shown
- Control soils have PFAS
- Conc increase with time for PFBS and PFOA in all trmts
- Superfund screening levels
  - PFBS 1600 mg/kg dry soil
  - PFOA 1.26 mg/kg dry soil
  - PFOS 1.26 mg/kg dry soil
- Some samples did not meet QA acceptance criteria
  - Biosolids controls 56 %
  - Solids application 23%
  - Control soil 8 %




# Other Observed PFAS

- Conc above RL are shown
  - Solids and control 30 ng/kg
  - Biosolids vary
- Biosolids show increasing conc with time
  - PFPeA
  - PFHxA
  - PFOA
- Solids show increasing conc with time
  - PFPeA
  - PFOA
- Control
  - Similar levels over time
  - Often similar to solids






# Oxidative Transformation to form PFOA



# PFAS Transformation Products


- Precursor conc range similar to PFAA conc
- Precursor conc decrease with time
- Stable PFAAs increased
- 85% mole balance

## Biosolids Controls



## Next Steps

- Complete data interpretation for study II
- Future Research
  - Characterize fate of precursors and transformation products
  - Characterize transport into vadose zone paired with ground water concentrations
  - Plan to use field and column studies





# Conclusions

- Metals
  - Sodium at background levels in 120 days
  - Copper conc in solids > control and liquids throughout the study
- NP
  - Liquids - removed after 120 days
  - Solids
    - Consistent with previous study, little change in conc for 1<sup>st</sup> 100 days
    - Slow decline throughout the study
  - Biosolids conc similar throughout the study
- PFAS
  - Observed in all samples
  - Lower Molecular Weight (MW) conc > higher MW conc
  - Precursors present and appear to convert to stable end products



# Acknowledgements

- Region 5 CRL: Francis Awanya, Anna Knoebel, Nidia Fuentes, Robert Thompson, Rob Snyder, Sylvia Griffin, Dennis Wesolowski, George Schupp
- PTSI: Solidea Bonina, Raghu Venkatapathy, Keith Bisbe, and Karen Koran
- ORD, NRMRL: Laura Boczek, Eric Rhodes, Marc Mills, Robert Ford, Jim Voit, Patrick Clark, Jonathan Ricketts (resigned EPA), Gerry Henderson, Michael Moeykens, Mark Kemper, Jennifer Goetz, Elizabeth Martin (resigned EPA), and Sean Murphy (resigned EPA)
- Others: Kavitha Dasu (currently Battelle Memorial Institute), Lawrence Wong (Senior Environmental Employee, CRL), Katrin Friesen (University of Alabama)



# Acronyms

- PFAS- per- and polyfluorinated alkyl substances
- PFCAs- perfluorinated carboxylic acids
- PFSAs- perfluorinated sulfonic acids
- PFHxA- perfluorohexanoic acid
- PFOA- perfluorooctanoic acid (MPFOA- isotopic version)
- PFOS- perfluorooctane sulfonic acid (MPFOS- isotopic version)
- PFHpA- perfluoroheptanoic acid
- PFPeA- perfluoropentanoic acid
- PFBS- perfluorobutane sulfonic acid
- PFHpS- perfluoroheptane sulfonic acid
- FTUCA- fluorotelomer unsaturated acid (8:2 measured)
- FTCA- fluorotelomer saturated acid (6:2, 8:2, 7:3 and 10:2 measured)
- RL – Reporting Limit
- WWTP- wastewater treatment plant
- MRM- multiple reaction monitoring
- RSD- relative standard deviation
- PFNA- perfluorononanoic acid (MPFNA- isotopic version)
- QA- quality assurance
- PFDA- perfluorodecanoic acid
- PFDoDA- perfluorododecanoic acid
- LC/MS/MS- liquid mass spectrometry
- GC/MS/MS- gas mass spectrometry
- PFAA- perfluorinated alkyl acid
- FTOH- fluorotelomer alcohol
- POTW- publicly owned treatment works
- MW- molecular weight



# Treatment Plots - November



**Liquid**




**Solid**




**Control**

# Field Study I NP results




# PFAS – More than just PFOA and PFOS




# PFAS analytes

## Perfluoroalkyl Carboxylates



PFBA	n = 4
PFPeA	n = 5
PFHxA	n = 6
PFHpA	n = 7
PFOA	n = 8
PFNA	n = 9
PFDA	n = 10
PFUdA	n = 11
PFDoA	n = 12
PFTTrA	n = 13
PFTeA	n = 14

## Perfluoroalkyl Sulfonates




PFBS	n = 4
PFPeS *	n = 5
PFHxS	n = 6
PFHpS	n = 7
PFOS	n = 8
PFNS *	n = 9
PFDS	n = 10




1611048\_22w Smooth(Mn,1x2)  
1611048-22

F14:MRM of 3 channels,ES-  
498.9 > 79.8  
1.502e+005



1611048\_22w Smooth(Mn,1x2)  
1611048-22

F14:MRM of 3 channels,ES-  
498.9 > 98.8  
1.239e+005



# PFAS Sampling

- **PFAS found in many common lab and field supplies and equipment**
  - Teflon - equipment, seals, sample caps, and bottles
  - Water proof paper and PPE
  - Personal care products
  - Clothing – water and stain repellent fabrics
  - Surface treatment on aluminum foil, food wrappers
  - Blue Ice
  - Supplies – sharpies, post-it notes
- **Avoid using these items when possible and pre-screen supplies and equipment**
  - Claims of PFOS/PFOA free may contain C6 versions
  - Read labels and product descriptions carefully
- **Information is evolving – check for updates**
- **Be careful about reusing existing equipment because of cross contamination – Decon and check for contamination**



# Equipment and supplies

## Avoid:

- Teflon, PTFE, and Fluoropolymers
- Aluminum foil may have PFAS surface treatment
- Decon 90, sharpies, post-it notes, waterproof papers or books
- Blue Ice
- Coated tyvek

## Acceptable

- HDPE, polypropylene, and silicone materials
- Alconox or Liquinox
- Ball point pens
- Water ice – double bag in polyethylene bags
- Uncoated Tyvek (if necessary)
- Sample bottles follow analytical SOP (usually PP or HDPE, not glass)



## Other precautions

- Food packaging may contain PFAS treatments – careful where you eat and wash hands before returning
- Frequent nitrile glove changes
- Collect sample, field, and equipment blanks
- Spiked blanks used by some



### Best practice

- Pretest materials and products for PFAS contamination
- Keep separate from “normal” supplies
- Test periodically for cross contamination